

# ELECTRONIC STETHOSCOPE



### Wireless Connection

Facilitates to remote auscultation and remote patient monitoring



### Improved sound quality

Enhanced sound volume and quality to improve auscultation




### Remotely Record and Share

With the help of voice recording app on your smart devices to save recordings for a second opinion



### Heart Sound and Lung Sound Filters

3 difference filter modes: heart sound, lung sound and wide range



Child



Adult



Type C



Audio jack 3.5mm

### Specification

Frequency response modes	Heart sound / Lung sound / Full frequency
Frequency range	Heart sound (20-200Hz) / Lung sound (100-1000Hz) / Full frequency (20-1200Hz)
Data wireless transmission	Bluetooth
Auscultation head type	Adult and Child sizes interchangeable with screw system. Compatible with 3M Littmann membrane.
Exterior dimensions	139 x 48 x 34 mm
Weight	190g
IP rating	IPX2
Power	2000 mAh Li-Polymer rechargeable battery
Battery life	1. 30 hours in auscultation mode 2. 200 hours in standby mode
Battery life recovery	Fully charged in 5 hours
Volume level	8 levels up to 50x amplification
Power consumption	90 mA

Bluetooth Application Scenario 1

Pairing with smart devices

Pairs with Smart Devices (Bluetooth 4.0 or later) for Auscultation Recording



Bluetooth 4.0  
or later



smartphones



tablets



laptops

1

Power on  
ES-2020 and  
Enter Bluetooth  
Pairing Mode



2

Turn on Bluetooth  
and select  
"ES-2020"



3

Open Recording  
App to Receive  
Signal from  
ES-2020




Bluetooth Application Scenario 2

Wireless monitoring

(Bluetooth earphones / speakers / personal sound amplification device)



Bluetooth  
transmitter



Bluetooth  
Speaker



Bluetooth  
earphones



personal sound  
amplification device

1

Connect Bluetooth  
Transmitter to 3.5  
mm Audio Jack



2


Turn on the Blue-  
tooth transmitter



Ensure Bluetooth Version  
Compatibility Before Use

3

Pair Bluetooth Audio  
Device according to  
instructions



Auscultation sound will be played through bluetooth  
device after successful pairing