

TD-1261C/D

EAR INFRARED THERMOMETER

FEATURES

- Fast ear scan in one second
- Dual-color backlight
- Probe cover ejection and detection
- Memory capacity up to 10 results
- Bluetooth data transmission function (TD-1261C)



THERMOMETER SPECIFICATIONS

Temperature Measurement Sites	Ear Temperature
Temperature Range	32°C ~ 43°C
Accuracy	±0.2°C (±0.4°F) for the measurement at 36°C-39°C (96.8°F-102.2°F) ±0.3°C (±0.5°F) outside above temperature range
Display Resolution	0.1°C or 0.1°F
Temperature Unit	°C or °F set by factory
Backlight Screen	Equipped (Normal temperature: Green backlight / Fever: Red backlight)
Fever/over Range Warning	Equipped
Probe Cover Detector	Equipped
Memory Capacity	10
Power Saving	3 min
Power Source	2 x AAA
Battery Life	1000 measurements
Dimension	160.5 (W) x 45.9 (H) x 32.5 (D) mm
Data Transmission	Bluetooth (TD-1261C); No Data Transmission (TD-1261D)
Weight	100 g (with battery)
Units per Carton	20
Net Weight per Carton	9.6 kg
Gross Weight per Carton	10.6 kg
Carton Size	46.7 (L) x 36.7 (W) x 24.2 (H) cm

PROBE COVER SPECIFICATIONS

Probe Cover Tray Dimension	148.76 (L) x 46 (W) x 21.7 (H) mm
Probe Cover Unit per Box	320 pcs
Probe Cover Box per Carton	6 pcs
Probe Cover Carton Size	31.5 (L) x 30 (W) x 15 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

6F., No.127, Wugong 2nd Rd., Wugong Dist.,
New Taipei City 24888, Taiwan
Tel: +886-2-6625-8188 Fax: +886-2-6625-0288