

# TD-4255

## BLOOD GLUCOSE MONITORING SYSTEM

7 Sec, 0.7 $\mu$ L

Strip-Ejection



### FEATURES

- LCD with backlight
- Strip-ejection function
- 7, 14, 21, 28, 60, 90 day average function



### SPECIFICATIONS

Enzyme Type	GDH-FAD
Sample Size	0.7 $\mu$ L
Reaction Time	7 seconds
Measurement Range	20 ~ 600mg/dL (1.1 ~ 33.3mmol/L)
Hematocrit Range	35-60%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100mg/dL (5.55 mmol/L); CV < 5% at 100mg/dL (5.55 mmol/L)
Accuracy	$\pm$ 15mg/dL if < 100mg/dL : $\pm$ 15% if $\geq$ 100mg/dL
Ketone Warning	Yes
Communication	Bluetooth(TD-4255C) / RS232(TD-4255D)
Power Source	AAA * 2
Memory Capacity	450 sets
Day Average	7, 14, 21, 28, 60, 90 days
Daily Alarm	4 Daily Alarms
Dimension	92.7 (L) x 60.8 (W) x 23.8 (H) mm
Weight	76.5g (Without Battery)
Operating Condition	10°C(50°F)~40°C(104°F)
Storage Condition	-20°C(-4°F) ~ +60°C(140°F) (Meter); 2°C(35.6°F)~32°C(89.6°F) (Strip)
Units Per Carton	30
Carton Size	405 (L) x 390 (W) x 365 (H) mm
Cross Weight Per Carton	7.5 kg



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,  
New Taipei City 24888, Taiwan

Tel : +886-2-6625-8188 Fax : +886-2-6625-0288