

TD-4287

BLOOD GLUCOSE MONITORING SYSTEM

FEATURES

- Compact size
- Before / After meal (AC/PC) in memory
- 150 memory storage



SPECIFICATIONS

Enzyme Type	GDH-FAD
Sample Size	0.5 μ L
Reaction Time	5 seconds
Measurement Range	20 ~ 600 mg/dL (1.1 ~ 33.3 mmol/L)
Hematocrit Range	35 - 60%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at 100 mg/dL (5.55 mmol/L)
Accuracy	\pm 15 mg/dL if < 100 mg/dL; \pm 15% if \geq 100 mg/dL
Communication	Mini USB
Power Source	CR2032 x 1
Memory Capacity	150 sets
Dimension	103.4 (L) x 28.1 (W) x 15.1 (H) mm
Weight	31.8 g (without battery)
Operating Condition	10°C (50°F) ~ 40°C (104°F)
Storage Condition	-20°C (-4°F) ~ +60°C (140°F) (Meter) 2°C (35.6°F) ~ 32°C (89.6°F) (Strip)
Units per Carton	40
Carton Size	520 (L) x 270 (W) x 363 (H) mm
Gross Weight per Carton	6 kg



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288