

# Vet Ca<sup>2+</sup>

TD-5220

## Plasma Calcium Monitoring System



### FEATURES

- Rapid Serum Calcium Test
- Quick Plasma Separation

### ACCESSORIES

Part no. 5220-CS001



Calcium test strip

Part no. 5220-BF002



Blood syringe filter

### SPECIFICATIONS

Sample Size	45 µL ±5 µL
Reaction Time	180 seconds
Measuring Range	0.5 to 3.0 mmol/L
Precision	≤ 1 mmol/L : SD < 0.12 mmol/L > 1 mmol/L : CV < 7.5%
Accuracy	±25%
Communication	RS232
Power Source	AAA x 2
Memory Capacity	300 sets
Dimension	166.1 (L) x 34.9 (W) x 34.9 (H) mm
Weight	81 g (without battery)
Operating Condition	2°C ~ 35°C
Storage Condition	-20°C (-4°F) ~ +60°C (140°F) (Meter); 2°C (35.6°F) ~ 30°C (86°F) (Strip)
Units per Carton	50
Carton Size	38 (L) x 21 (W) x 35.5 (H) cm
Gross Weight per Carton	10.5 kg



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,  
New Taipei City 24888, Taiwan  
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288