

TD-7014

PORTABLE COMPRESSOR NEBULIZER

DUAL POWER OPTIONS



Battery operated



Power adaptor

DC MICROPUMP TECHNOLOGY



Low-noise




Handheld

EFFICIENT

Fine particle size and effective dose delivery

Dual filter mechanism to filter out the small particles

The ring structure shortens the nebulization process and enhances the nebulization efficiency



SPECIFICATIONS

Particle Size	Approx. 3 microns
Nebulization Rate	≥ 0.25mL/min
Operating Pressure Range	35~69kPa (5~10 Psi)
Medication Capacity	6mL
Power Consumption	6W
Power supply	AC adaptor or 4 x AA battery (Alkaline or rechargeable)
Waterproof grade	IP22
Dimension	115L / 85W / 207H (mm)
Weight	205g (without batteries)
Units Per Carton	24 PCS
Carton Size	L523 x W493 x H376 mm
Gross Weight Per Carton	12 KGS



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288