

TD-3124

BLOOD PRESSURE MONITORING SYSTEM



FEATURES

- 100 sets of memory
- IHB detection
- One button and large LCD display design; simple operation

SPECIFICATIONS

Accuracy of Pressure	± 3 mmHg or $\pm 2\%$ of reading
Accuracy of Heart Rate	$\pm 4\%$ of reading
Measurement Range	Systolic : 60 ~ 255 mmHg / Diastolic : 30 ~ 195 mmHg Pulse rate : 40~199 / minute
Power Source	4 x AA or 6V DC adaptor (optional)
Dimension	130.5 (L) x 90 (W) x 55 (H) mm
Weight	210 g (without battery)
Operating Condition	5 °C ~ 40 °C ; Humidity 15% ~ 93%
Battery Life	At least 500 measurements
Power Saving	180 seconds
Units per Carton	30
Net Weight per Carton	10.8 kg
Gross Weight per Carton	11.8 kg
Carton Size	42.5 x 32.5 x 34.5 cm

ACCESSORY

Part. No.	Cuff Size	Arm Circumference
CF3312	Medium	24~35cm
CF3314	Wide Range	24~43cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288