

TD-3128

BLOOD PRESSURE MONITORING SYSTEM







Data Transmission

4 Users Selection



FEATURES

-  Average 3 measurements
-  Large screen 121 x 63mm
-  IHB Detection
-  Low-noise



Healthy Check



Healthy Check Pro

SPECIFICATIONS

Accuracy of Pressure	± 3 mmHg or $\pm 2\%$ of reading
Accuracy of Heart Rate	$\pm 4\%$ of reading
Measurement Range	Systolic : 60 ~ 255 mmHg / Diastolic : 30 ~ 195 mmHg Pulse rate : 40~199 / minute
Memory Capacity	400 sets
Data Output	Bluetooth (TD-3128B) / RS232(TD-3128C)
User Number	4 sets
Power Source	4 x AA or 6V DC Adaptor (Optional)
Dimension	141 (L) x 121 (W) x 72 (H) mm
Weight	350 g (without battery)
Operating Condition	5 °C ~ 40 °C ; Humidity 15% ~ 93%
Battery Life	At least 200 measurements
Power Saving	180 seconds
Units per Carton	12
Net Weight per Carton	15.7 kg
Gross Weight per Carton	16.7 kg
Carton Size	47.3 x 35.8 x 41.5 cm

ACCESSORY

Part No.	Cuff Size	Arm Circumference
CF3312	Medium	24-35cm
CF3314	Wide Range	24-43cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288