


Certificate

Taidoc Technology Corporation

using

Immunoassay on TD-4611

participated in the IFCC HbA1c Certification Programme to demonstrate traceability to the IFCC Reference Measurement Procedure and performed as shown below.



Total Error =	2.5	mmol/mol
Bias =	-0.1	mmol/mol
Imprecision =	2.4	%
Grade =	Bronze	

Criteria derived from the IFCC model for Quality Targets HbA1c (Clin Chem 2015;61 : 752-59)

Date of Certification : 01 January 2020

Date of Expiry : 01 January 2021


 IFCC Network Coordinator
 Dr. C.W. Weykamp