


TD-4140B

LIGHTNING DONGLE



Blood Glucose
Ketone
Cholesterol



Compatible Device

- iPod touch 6
- iPad Air, iPad Air 2
- iPad (5th/6th generation)
- iPhone 5s, iPhone 6 and later
- iPad mini 2, iPad mini 3, iPad mini 4



SPECIFICATIONS / Strip

	Strip I	Strip II	Strip III
Parameter	Glucose	Ketone	Cholesterol
Sample Size	1.0 μ L	1.0 μ L	3.0 μ L
Reaction Time	5 seconds	10 seconds	\leq 60 seconds
Measurement Range	10 ~ 700 mg/dL (0.56 ~ 38.89mmol/L)	0.1 ~ 8.0 mmol/L	100 ~ 400 mg/dL (2.6 ~ 10.4 mmol/L)
Hematocrit Range	0 ~ 70 %	10 ~ 70%	20 ~ 60%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at 100 mg/dL (5.55 mmol/L)	\leq 1 mmol/L, SD < 0.1 mM; > 1 mmol/L, CV<7.5%	CV < 7.5%
Package	Vial pack	Single foil pack	Single foil pack

SPECIFICATIONS / Lightning Dongle

Communication	Lightning
Dimension	51 (L) \times 24 (W) \times 8.8 (H) mm
Weight	9.0 g
Operating Condition	+10°C (50°F) ~ +40°C (104°F)
Storage Condition	-20°C (-4°F) ~ +60°C (140°F)
Units per Carton	90 pcs
Carton Size	43 (L) \times 35 (W) \times 39.5 (H) cm
Gross Weight per Carton	12.6 kg



TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288

www.taidoc.com