

TD-7014 JET NEBULIZER



Handheld



Low-noise



Battery operated / Power adaptor / Power bank cable



SPECIFICATIONS



Particle Size MMAD	Approx. 2.5 μ m
Nebulization Rate	\geq 0.25mL/min
Operating Pressure Range	35-69kPa (5-10 Psi)
Medication Capacity	10mL
Power Consumption	6W
Power supply	AC adaptor or 4 x AA battery (Alkaline or rechargeable)
Compressor Pressure Range	90 to 138 kPa
Flow Rate	\geq 3 LPM(L/min)
Operating Cycle	30min on / 30min off
Waterproof grade	IP22
Dimension	88L / 85W / 238H (mm)
Weight	195g (without batteries)
Units Per Carton	24 PCS
Carton Size	L523 x W493 x H376 mm
Gross Weight Per Carton	12 KGS




TaiDoc Technology Corp.





B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel: +886-2-6625-8188 Fax: +886-2-6625-0288

www.taidoc.com

Medication Cup 10
MdC 10



Accessories

Part no.	NEB-ACC-001	NEB-ACC-002	NEB-ACC-003	NEB-ACC-004
Description	MdC 10 kit	MdC 11 kit	Adult Mask	Child Mask
Picture				



TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288

www.taidoc.com