

# TD-8255 Fingertip Pulse Oximeter



Silicone Rubber


 Bluetooth®



TD-8255A



TD-8255B



Large LCD with backlight



Perfusion Index

## SPECIFICATIONS

Part No.	TD-8255A(No Data Transmission) / TD-8255B(Bluetooth)
Communication	Bluetooth
Dimension	63 (L) x 35 (W) x 32 (H) mm
Weight	40g (without battery)
Heart Rate Range	30 ~ 250 beats per minutes (bpm)
SpO2 Range	0% ~ 100%
Accuracy (SpO2)	80%~100% ±2 digits
Method	Dual wavelength LED
Perfusion Index	5-bar indicator
Accuracy (Heart rate)	±1 bpm or ±1% whichever is greater
Power Source	2 x AAA
Time to power saving	15 seconds after removing finger
Backlight	LCD backlight
Operating Temperature	10°C ~ 40°C
EMC Standard	IEC60601-1-2
Safety Standard	IEC60601-1
Performance Standard	ISO80601-2-61
Unit per Carton	80
Net Weight per Carton	11.2kgs
Gross Weight per Carton	12.5kgs
Carton Size	518 (L) x 423 (W) x 300 (H) mm



iOS

Android



Healthy Check



iOS

Android



Healthy Check Pro



TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,  
New Taipei City 24888, Taiwan  
Tel: +886-2-6625-8188 Fax: +886-2-6625-0288

www.taidoc.com