


TD-3261

Blood Pressure and Blood Glucose Monitor



Telemedicine without time difference!



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No. 127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288

Blood Pressure and Blood Glucose Monitor **2in1**



Blood Pressure monitor



Clever inflation technology



Irregular Heartbeat (IHB) detection



Automated Averaging Function



Smart Auto-Cycling Function



Blood Glucose monitor



No-coding blood glucose test strip



Pre-meal (AC) and Post-meal (PC) records



Strip ejection



SPECIFICATIONS

Power Source	Rechargeable Li-Ion Battery and External DC 6V
Memory Capacity	Maximum 864 memory records
Communication	Bluetooth smart (TD-3261B) / 4G (TD-3261G)
Power Saving	180 seconds
Size of Monitor w/o Cuff	147mm (L) x 105mm (W) x 80mm (H), 500g (include battery)
Operating Condition	10°C~40°C

BLOOD GLUCOSE MONITOR

Enzyme type	GOD
Sample Size / Reaction Time	0.5 μ L / 5 seconds
Measurement Range	20-600 mg/dL (1.1-33.3 mmol/L)
Hematocrit Range	20% ~ 60%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100mg/dL (5.55 mmol/L); CV < 5% at 100mg/dL (5.55 mmol/L)
Accuracy	\pm 15mg/dL if \leq 100 mg/dL; \pm 15% if > 100 mg/dL
Ketone Warning	Yes

BLOOD PRESSURE MONITOR

Accuracy of Pressure	\pm 3 mmHg or \pm 2% of reading
Accuracy of Heart Rate	\pm 4% of reading
Measurement Range	Systolic : 60 ~ 250 mmHg / Diastolic : 30 ~ 180 mmHg Pulse rate : 40 ~ 199 / minute
Units per Carton	12
Net Weight per Carton	11.36 kg
Gross Weight per Carton	12.36 kg
Carton Size	L42.5 x W37.0 x H34.5 cm

ACCESSORY

Part. No.	Cuff Size	Arm Circumference
CF3312	Meium	24-35 cm
CF3314	Wide Range	24-43 cm