

TD-4255

BLOOD GLUCOSE MONITORING SYSTEM

7 Sec, 0.7 μ L

Strip-Ejection



FEATURES

- LCD with backlight
- Strip-ejection function
- 7, 14, 21, 28, 60, 90 day average function



SPECIFICATIONS

Enzyme Type	GDH - FAD
Sample Size	0.7 μ L
Reaction Time	7 seconds
Measurement Range	20 ~ 600 mg/dL (1.1 ~ 33.3 mmol/L)
Hematocrit Range	35 - 60%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at 100 mg/dL (5.55 mmol/L)
Accuracy	\pm 15 mg/dL if < 100 mg/dL; \pm 15% if \geq 100 mg/dL
Ketone Warning	Yes
Communication	Bluetooth (TD-4255C) / RS232 (TD-4255D)
Power Source	AAA * 2
Memory Capacity	450 sets
Day Average	7, 14, 21, 28, 60, 90 days
Daily Alarm	4 Daily Alarms
Dimension	92.7 (L) x 60.8 (W) x 23.8 (H) mm
Weight	76.5g (Without Battery)
Operating Condition	10°C (50°F) ~ 40°C (104°F)
Storage Condition	-20°C (-4°F) ~ +60°C (140°F) (Meter); 2°C (35.6°F) ~ 32°C (89.6°F) (Strip)
Units Per Carton	30
Carton Size	405 (L) x 390 (W) x 365 (H) mm
Gross Weight Per Carton	7.5 kg



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288