

TD-1090

Thermometer Chamber



FEATURES

- Accuracy verification for ear and forehead thermometers
- Adjustable target temperature
- Data export capability

Thermometer Fixter

Calibration Cable

Part no.	ACC-1090-001	ACC-1090-002	ACC-1090-003	ACC-1090-005	ACC-1090-004
Picture	  	 	 	 	 
	TD-1107 TD-1261	TD-1241	TD-1115	TD-1242	TD-1107

SPECIFICATIONS

Display Resolution	0.1°C/0.1°F
Power & Fuse(s)	100 ~ 240VAC (50/60Hz) / 3A
Memory	1000
Operating Condition	16°C to 32°C (60.8°F to 89.6°F), 85% RH or less
Temperature Range (0.98 emissivity)	37.0°C ~ 41.0°C (98.6°F ~ 105.8°F)
Acceptable Recalibration Range	±0.5°C (±0.9°F)
Stability	±0.03°C at 37°C
Measurement Units	°C (default) or °F
Warm-up Time at 25°C Ambient	30 minutes
Nominal Emissivity	0.98 emissivity
Computer Interface	USB HID
Dimension	300mm x 350mm x 160mm

SHIPPING INFO

Units per Carton	1
Gross Weight per Carton	6.93 kg
Carton size	41.3 (L) x 22.3 (W) x 37.6 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel: +886-2-6625-8188 Fax: +886-2-6625-0288