

TD-1115

Non-Contact Forehead Thermometer



On/Scan Button

Display Screen

Forehead/
Surface Mode



FEATURES

- Non-Contact infrared Thermometer
- Forehead and Surface temperature
- 3 Colors Backlight Indicator

SPECIFICATIONS

Measurement Mode	Forehead	Surface
Measurement Range	34.4°C ~ 42.2°C (93.9°F ~ 108°F)	1°C ~ 60°C (33.8°F ~ 140°F)
Backlight indicator	34.4°C 38°C 42.2°C Lo Normal Fever Hi	1°C 60°C Lo Range Hi
Accuracy	35°C ~ 42°C (± 0.2°C) < 35°C or > 42°C (± 0.3°C)	35°C ~ 42°C (± 0.2°C) < 35°C or > 42°C (± 0.3°C)
Operating Condition	10°C ~ 40°C (50°F ~ 104°F) 30-85% R.H. (non-condensing)	
Storage/Transportation Condition	-20°C ~ 60°C (-4°F ~ 140°F) 30-85% RH. (non-condensing)	
Resolution	0.1°C/0.1°F	
Measurement Unit	°C or °F	
Power Source	CR2032 x 1	
Battery Life	At least 1000 measurements	
Power Saving	30 seconds	
Dimension	127.7(L) x 37.7(W) x 20.7(H) mm	
Weight	56.2 g	

SHIPPING INFO

Units per Carton	80
Gross Weight per Carton	11.2 kg
Carton size	43.6 (L) x 26.5 (W) x 42 (H) cm



TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288

www.taidoc.com