

TD-1241

Non-Contact Forehead Thermometer



FEATURES

- 3 Measuring mode:
Forehead/Surface/Room temperature
- Accurate and Fast
- LCD Backlight



Healthy Check



Healthy Check Pro

SPECIFICATIONS

Dimension	150 (L) x 48.5 (W) x 55 (H)mm
Weight	125.8 g (With Battery)
Resolution	0.1°C/0.1°F
Measurement Range	Forehead: 22°C-44°C Surface: 0°C -100°C Room: 10°C -40°C
Accuracy	Forehead: 36°C-39°C±0.2°C(22°C-35.9°C, 39.1°C -44°C±0.3°C) Surface/Room: ±.1°C
Communication	No Data Transmission (TD-1241A) / Bluetooth 4.0 (TD-1241B)
Measurement Units	°C(preset)/ °F
Power source	2*1.5V AA
Battery life	> 3000 measurements
Memory	30 measurements
Power saving	30 seconds
Operating condition	10°C to 40°C (50°F to 104°F), 95% RH or less (non-condensing)
Storage/Transportation condition	-20°C to 60°C (-4°F to 140°F), 95% RH or less (non-condensing)
Measuring Distance	3-7cm

SHIPPING INFO

Units per Carton	30
Gross Weight per Carton	13.5 kg
Carton Size	42 (L) x 39 (W) x 41 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel: +886-2-6625-8188 Fax: +886-2-6625-0288