

# TD-2555

## WEIGHT SCALE

**Mb<sup>H</sup>**  
Mobile Health



 **Bluetooth®**



Body mass index function



Age, height, and gender for more precise algorithm



Large LCD with backlight



Automatically power-off function



iOS

Android



Healthy Check



iOS

Android



Healthy Check Pro

### SPECIFICATIONS

BMI (Body Mass Index)	Yes
Max. capacity	250 Kg
Accuracy	4.0 kg ~ 60.0 kg ± 0.3 kg, 60.1 kg ~ 250.0 ± 0.5%
Graduation	0.1 Kg
Age Range	10 to 80
Height Range	100 to 220 cm
User	5
Communication	Bluetooth smart
Memory	135 sets
Backlight	LCD Backlight
Power Source	4 x AAA
Dimension	340(L) x 400 (W) x 26 (H) mm

### SHIPPING INFO

Units per Carton	5
Gross Weight per Carton	18 kg
Carton Size	49.3 (L) x 29.3 (W) x 45.5 (H) cm



www.taidoc.com

**TaiDoc Technology Corp.**

B1-7F, No.127, Wugong 2nd Rd., Wugu Dist.,  
New Taipei City 24888, Taiwan  
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288