

TD-3127

BLOOD PRESSURE MONITORING SYSTEM

FEATURES

- Arm type
- IHB detection
- 4 users selection
- Large LCD display

4 Users Selection



SPECIFICATIONS

Accuracy of Pressure	± 3 mmHg or $\pm 2\%$ of reading
Accuracy of Heart Rate	$\pm 4\%$ of reading
Measurement Range	Systolic : 60 ~ 255 mmHg / Diastolic : 30 ~ 195 mmHg Pulse rate : 40~199 / minute
Memory Capacity	100 sets
User Number	4 sets
Power Source	4 x AA or 6V DC Adaptor (Optional)
Dimension	150 (L) x 96 (W) x 66 (H) mm
Weight	150 g (without battery)
Operating Condition	5 °C ~ 40 °C ; Humidity 15% ~ 93%
Battery Life	At least 300 measurements
Power Saving	180 seconds

SHIPPING INFO

Units per Carton	20
Gross Weight per Carton	12 kg
Carton Size	56.8 (L) x 36.8 (W) x 28.7 (H) cm

ACCESSORY

Part. No.	Cuff Size	Arm Circumference
CF3312	Medium	24~35cm
CF3314	Wide Range	24~43cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288