

TD-3140

UPPER ARM BLOOD PRESSURE MONITOR

Mb[®]H
Mobile Health



Smartphone App for cardiovascular evaluation

-  **HI** Fast Heart Beats
Heart rate: > 110 bpm.

-  **LO** Slow Heart Beats
Heart rate: < 50 bpm.

-  **IHB** Irregular Heart Beats
Heart rate: premature, delayed contraction, skipped, pause between heartbeats, extra systole.

-  **AF** Atrial Fibrillation
Heart rate: > 100 bpm, quivering or irregular heartbeat.

-  **Cardiovascular Evaluation**
Evaluation: cardiovascular or age assessment.





MbH BP



SPECIFICATIONS

Accuracy of Pressure	± 3 mmHg or $\pm 2\%$ of reading
Accuracy of Heart Rate	$\pm 4\%$ of reading
Measurement Range	Systolic: 60 ~ 255 mmHg / Diastolic: 30 ~ 195 mmHg Pulse rate: 40~199 / minute
Memory Capacity	200 sets
Communication	Bluetooth 4.0
Power Source	4 x AAA
Dimension	145 (L) x 64.5 (W) x 29.6 (H) mm
Weight	240 g (without battery)
Operating Condition	5 °C ~ 40 °C; Humidity 15% ~ 93%
Battery Life	At least 200 measurements
Power Saving	180 seconds
Backlight	LCD backlight
Net Weight per Carton	12.63 kg
Cuff Size	Wide Range: 24 ~ 43 cm, M Size: 24 ~ 35 cm

SHIPPING INFO

Units per Carton	30
Gross Weight per Carton	15.6 kg
Carton Size	42 (L) x 39 (W) x 41 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F, No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288