

TD-4121

BLOOD GLUCOSE MONITORING SYSTEM



±15mg/dL & ±15%

5 Sec, 0.5 µL

Touchscreen

Strip-Ejection



FEATURES

- Strip port indicator
- Dual language displays, English and Spanish
- Before / After meal (AC/PC) presetable
- Emergency call button



SPECIFICATIONS

Enzyme Type	GDH-FAD
Sample Size	0.5 µL
Reaction Time	5 seconds
Measurement Range	20 ~ 600 mg/dL (1.1 ~ 33.3 mmol/L)
Hematocrit Range	35% - 60%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at 100 mg/dL (5.55 mmol/L)
Accuracy	±15mg/dL if < 100 mg/dL; ±15% if ≥ 100 mg/dL
Ketone Warning	Yes
Communication	2.5G (GPRS, EDGE) + 3G (UMTS, HSPA, HSPA+)
Power Source	Li-Ion rechargeable
Memory Capacity	450 sets
Day Average	7, 14, 21, 28, 60, 90 days
Daily Alarm	4 daily alarms
Dimension	99 (L) x 59 (W) x 15.6 (H)mm
Weight	51.4 g (without battery)
Operating Condition	10°C (50°F) ~ 40°C (104°F), below 85% R.H.
Storage Condition	-20° C (-4°F) ~ 60°C (140°F) (Meter) 2°C (35.6°F) ~ 30°C (86.0°F) (Strip)
Units per Carton	40
Carton Size	46.7 (L) x 29.4 (W) x 41.5 (H) cm
Gross Weight per Carton	11.84 kg



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288