

TD-4235A

BLOOD GLUCOSE MONITORING SYSTEM



±15mg/dL & ±15%

USB Port

7 Sec, 0.7 µL



FEATURES

- No coding required
- Before / After meal (AC/PC) presettable
- LCD with backlight
- Alternative Site Testing (AST)
- Display calendar and clock when meter standby

SPECIFICATIONS

Enzyme Type	GOD
Sample Size	0.7 µL
Reaction Time	7 seconds
Measurement Range	20 ~ 600 mg/dL (1.1 ~ 33.3 mmol/L)
Hematocrit Range	30% - 55%
Precision	±5% with respect to Standard
Accuracy	±15 mg/dL if ≤ 100 mg/dL; ±15% if > 100 mg/dL
Ketone Warning	No
Communication	USB
Power Source	AAA x 1
Memory Capacity	150 sets
Day Average	No
Daily Alarm	4 daily alarms
Dimension	93 (L) x 36 (W) x 21.3 (H) mm
Weight	45.4 g (with battery)
Operating Condition	10°C ~ 40°C, below 85% R.H.
Storage Condition	-20 °C ~ +60 °C (Meter); 4°C ~ 40°C (Strip)



SHIPPING INFO

Units per Carton	40
Gross Weight per Carton	11.2 kg
Carton Size	52.8 (L) x 27.8 (W) x 38.3 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288