

TD-4258

POCT KETONE MONITORING SYSTEM



FEATURES

- Strip eject function
- USB & WIFI data transmission
- Embedded barcode scanner
- Flexible Memory Capacity: 1000 Measurement Results, 2000 Operator ID, 500 QC



SPECIFICATIONS

| | |
|---------------------|---|
| Enzyme Type | HBD |
| Sample Size | 1.0 μ L |
| Reaction Time | 10 seconds |
| Measurement Range | 0.1 ~ 8.0mmol/L (1.04 ~ 83.2mg/dL) |
| Hematocrit Range | 10 ~ 70% |
| Precision | ≤ 2 mmol/L, SD < 0.2mM : >2mmol/L, CV < 7.5% |
| Accuracy | ± 0.5 mmol/L if ≤ 2 mmol/L : $\pm 25\%$ if > 2mmol/L |
| Communication | WiFi & Mini-USB |
| Power Source | Li-ion Rechargeable Battery |
| Memory Capacity | 1000 Measurement Results, 2000 Operator ID, and 500 QC |
| Dimension | 130.47 (L) x 65.39 (W) x 35.48 (H) |
| Weight | 170.8 g with battery |
| Operating Condition | 10°C (50°F) ~ 40°C (104°F) |
| Storage Condition | -20°C (-4°F) ~ 60°C (140°F) (Meter); 2°C (35.6°F) ~ 32°C (89.6°F) (Strip) |

SHIPPING INFO

| | |
|-------------------------|-----------------------------------|
| Units Per Carton | 12 |
| Gross Weight Per Carton | 12 kg |
| Carton Size | 43.3 (L) x 37.8 (W) x 36.5 (H) cm |

Docking Station

Part no. 9019-DS02

- Wireless charging in compliance with WPC1.1
- Dimension 88.99 x 190.55 x 121.54 mm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288