

TD-4261A

LACTATE MONITORING SYSTEM



FEATURES

- Strip eject function
- Mini USB communication
- TFT color screen



mini USB communication



lightweight



compact size



Strip Ejection

0.8µL, 5 sec

SPECIFICATIONS

Sample Size	0.8 µL
Reaction Time	5 seconds
Measurement Range	0.3-22 mmol/L
Hematocrit Range	10% - 65%
Precision	≤3mmol/L, SD<0.3mM; >3mmol/L, CV<7.5%
Package	Single foil pack
Communication	Mini USB port
Power Source	Lithium-ion battery rechargeable
Memory Capacity	450 sets
Dimension	96 (L) x 61 (W) x 26 (H) mm
Weight	67.2 g
Operating Condition	strip / meter: +10°C (50°F) ~ 40°C (104°F)
Storage Condition	strip: 2°C (35.6°F) ~ 30°C (86°F) / meter: -20°C (-4°F) ~ 60°C (140°F)



Lactate

SHIPPING INFO

Units Per Carton	40
Gross Weight Per Carton	12.8 kg
Carton Size	49.8 (L) x 35.8 (W) x 36 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel: +886-2-6625-8188 Fax: +886-2-6625-0288