

# TD-4283

## Blood Glucose monitoring

### FEATURES

- Compact size
- Easy to use
- 150 memory sets



### SPECIFICATIONS

Enzyme type	GDH-FAD
Sample Size	0.5 µL
Reaction Time	5 seconds
Measurement Range	20 ~ 600 mg/dL (1.1~33.3 mmol/L)
Hematocrit Range	35% - 60%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at 100 mg/dL (5.55 mmol/L)
Accuracy	±15 mg/dL if < 100 mg/dL; ±15% if ≥ 100 mg/dL
Communication	Strip Port Communication
Power Source	CR2032 x 1
Memory Capacity	150 sets
Dimension	43.5 (L) x 39.3 (W) x 24.2 (H) mm
Weight	13 g (without battery)
Operating Condition	+10°C (50°F) ~ +40°C (104°F), below 85% R.H.
Storage Condition	-20°C (-4°F) ~ +60°C (140°F) (Meter) ; 2°C (35.6°F) ~ 32°C (89.6°F) (Strip)

### SHIPPING INFO

Units per Carton	70
Gross Weight per Carton	12.6 kg
Carton Size	59.3 (L) x 33.8 (W) x 31.7 (H) cm

### SPECIFICATIONS-Accessories

ACC-001	Bluetooth dongle
ACC-002	USB cable



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,  
New Taipei City 24888, Taiwan  
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288