


TD-7013

COMPRESSOR

NEBULIZER

Medication Cup Technology



FEATURES

- DC micropump technology
- Compact size for better usability
- Fine particle size for better treatment
- One step operation, simple and quiet



Low-noise

SPECIFICATIONS

Particle Size	MMAD approx. 3 um
Nebulizer Rate	>0.35 mL/min
Medication Capacity	6mL
Power Consumption	12W
Power Source	12V DC (Switch power universal adapter)
Dimension	125 x 115 x 46 mm
Weight	220g
Operating Condition	5°C ~ 40°C; RH 15% ~ 93%
Storage Condition	-25°C ~ 70°C; RH 10% ~ 95%
Net Weight per Carton	12 kg

SHIPPING INFO

Units per Carton	20
Gross Weight per Carton	16 kg
Carton Size	57.8 (L) x 33.8 (W) x 41 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288