



CHOLESTEROL

GLUCOSE

KETONE

**TRIPLE IN ONE**

Use only with Procheck Witty Blood Glucose /  $\beta$ -Ketone / Total Cholesterol Test Strips



**Pocket-Size**  
No need for a bulky case



**ProCheck APP**  
ProCheck app for personal health management



**Battery Free**  
The power comes from the Smartphone


## Procheck Witty Smartphone Dongle



**USB-C Connector**  
Connects directly to your Smartphone



**Lightning**  
Connects directly to your Smartphone



## Procheck Smart Smartphone Dongle



## SPECIFICATIONS / STRIP

### GLUCOSE



Sample Size	0.8 $\mu$ L
Reaction Time	5 seconds
Measurement Range	10 ~ 700 mg/dL (0.56 ~ 38.89 mmol/L)
Hematocrit Range	0 ~ 70 %
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at $\geq$ 100 mg/dL (5.55 mmol/L)
Package	Vial Pack (50 pcs x 1 vial)
Package Dimension	43 (L) x 40 (W) x 58 (H) mm
Shipping INFO	
Units per Carton	400
Gross Weight per Carton	15.9 kg
Carton Size	444 (L) x 393 (W) x 402 (H) mm

### KETONE



Sample Size	1.0 $\mu$ L
Reaction Time	10 seconds
Measurement Range	0.1 ~ 8.0 mmol/L
Hematocrit Range	10 ~ 70 %
Precision	$\leq$ 1 mmol/L, SD < 0.1 mM; > 1 mmol/L, CV < 7.5%
Package	Single Foil Pack (50 pcs)
Package Dimension	96 (L) x 10 (W) x 220 (H) mm
Shipping INFO	
Units per Carton	720
Gross Weight per Carton	NA
Carton Size	456 (L) x 378 (W) x 368 (H) mm

## Android Compatible Device (Recommends operating system 8.0 or later)

1. Samsung Galaxy A70 / A8+ / S8 / S9 / S9+ / S10 / S10+ / Note 8 / Note 9 / Note 10
2. Xiaomi Mi 9T / Mi 11 Lite / Redmi Note 7 / Redmi Note 8 Pro
3. Google Pixel 3 / Pixel 4 / Pixel 5
4. OPPO Reno Z 5. Realme XT 6. Nokia 8.1

\*Some android phones might not be compatible



Google Play and the Google Play logo are trademarks of Google LLC.

## iOS Compatible Device (Requires operating system 12.0 or later)

1. iPad Air, iPad mini 3
2. iPad (6th generation), iPod touch (6th generation)
3. iPhone 5s / 6 / 6s Plus / 7 Plus / XR / XS / 11 / 12 Pro Max

\*Some ios phones might not be compatible



Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries and regions. App Store is a service mark of Apple Inc.

## SPECIFICATIONS

Part No.	Procheck Witty
Communication	USB-C
Dimension	55 (L) X 40.29 (W) X 10.35 (H) mm
Weight	11.8 g
Operating Condition	+10°C (50°F) ~ +40°C (104°F)
Storage Condition	-20°C (-4°F) ~ +60°C (140°F)
Shipping INFO	
Units per Carton	30
Gross Weight per Carton	8.5 kg
Carton Size	340 (L) x 290 (W) x 383 (H) mm

Part No.	Procheck Smart
Communication	Lightning
Dimension	51 (L) X 24 (W) X 8.8 (H) mm
Weight	9 g
Operating Condition	+10°C (50°F) ~ +40°C (104°F)
Storage Condition	-20°C (-4°F) ~ +60°C (140°F)
Shipping INFO	
Units per Carton	30
Gross Weight per Carton	8.5 kg
Carton Size	340 (L) x 290 (W) x 383 (H) mm

## CHOLESTEROL



Sample Size	3.0 $\mu$ L
Reaction Time	$\leq$ 60 seconds
Measurement Range	100 ~ 400 mg/dL (2.6 ~ 10.4 mmol/L)
Hematocrit Range	20 ~ 60%
Precision	CV < 7.5%
Package	Single Foil Pack (10 pcs)
Package Dimension	96 (L) x 5.5 (W) x 220 (H) mm
Shipping INFO	
Units per Carton	280
Gross Weight per Carton	NA
Carton Size	456 (L) x 378 (W) x 368 (H) mm