



TD-4257

BLOOD GLUCOSE MONITORING SYSTEM

FEATURES

- EF (Electricity Fractionation) technology – HCT range **20-65%**
- Target range indicator



SPECIFICATIONS

Enzyme Type	GDH-FAD
Sample Size	0.5 μ L
Reaction Time	5 seconds
Measurement Range	20 ~ 650 mg/dL (1.1~36.1mmol/L)
Hematocrit Range	20% ~ 65%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at 100 mg/dL (5.55 mmol/L)
Accuracy	± 15 mg/dL if < 100 mg/dL; $\pm 15\%$ if ≥ 100 mg/dL
Communication	Strip port
Ketone Warning	Yes
Power Source	2 x AAA
Memory Capacity	400 sets
Day Average	7, 14, 21, 28, 60, 90 days for blood glucose
Daily Alarm	4 daily alarms
Dimension	97 (L) x 53.7 (W) x 21.7 (H)mm
Weight	81g (without battery)
Operating Condition	10°C (50°F) ~ 40°C (104°F), 10% ~ 85% R.H.
Storage Condition	Meter: -20°C (-4°F) ~ 60°C (140°F); Strip: 2°C (35.6°F) ~ 30°C (86°F)
Package	Vial pack / Single foil pack

SHIPPING INFO

Units per Carton	40
Gross Weight per Carton	12.8 kg
Carton size	49.8 (L) x 35.8 (W) x 36 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288