


TD-5083

AKTIVSAFE LANCET




FEATURES


- Easy to use
- Single-use only
- Lancet automatic retraction



1
Twist off the protective cap and remove the cap



2
Press the safety lancet body firmly against the puncture site to activate the safety lancet



3
Dispose the used safety lancet and collect blood sample



Lancet

Lancet Gauge	28G
Penetration Depth	1.5±0.2mm
Packaging	100pcs/Box
Units per Carton	6
Gross Weight per Carton	2.5kg
Carton Size	30.5 (L) x29 (W) x13.5 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel: +886-2-6625-8188 Fax: +886-2-6625-0288