


Certificate

Taidoc Technology Corporation

using

Immunoassay on TD-4611

participated in the IFCC HbA1c Certification Programme to demonstrate traceability to the IFCC Reference Measurement Procedure and performed as shown below.




Total Error =	2.6	mmol/mol
Bias =	+1.2	mmol/mol
Imprecision =	1.4	%
Grade =	Bronze	

Criteria derived from the IFCC model for Quality Targets HbA1c (Clin Chem 2015;61 : 752-59)

Date of Certification : 01 January 2022

Date of Expiry : 01 January 2023


IFCC Network Coordinator
C. Siebelder