


TD-5052

AKTIVSAFE LANCET




FEATURES


- Easy to use
- Five different size
- Single-use only
- Lancet automatic retraction
- FDA 510(k) cleared Class II medical device.



1 Twist off the protective cap and remove the cap



2 Press the safety lancet body firmly against the puncture site to activate the safety lancet



3 Dispose the used safety lancet and collect blood sample

Gauge	21G	23G	26G	28G	30G
Lancet color	 Green	 Blue	 Yellow	 Orange	 Purple
Penetration Depth	2.2 mm	2.2 mm	1.8 mm	1.8 mm	1.6 mm
Packaging	100 pcs/box 60 boxes/carton	100 pcs/box 60 boxes/carton	100 pcs/box 60 boxes/carton	100 pcs/box 60 boxes/carton	100 pcs/box 60 boxes/carton

SHIPPING INFO

Units per Carton	60
Gross Weight per Carton	18 kg
Carton Size	46.3 (L) x 38.7 (W) x 43.5 (H) cm



TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288

www.taidoc.com