

TD-1035

Wearable Bluetooth Thermometer

Mb^H
Mobile Health



SPECIFICATIONS

Measurement Range	20.00°C (68.00°F) ~ 45.00°C(113.00°F)
Accuracy	20.00°C ~ 34.99°C (68.00°F ~ 94.98°F): ±0.10°C(0.18°F) 35.00°C ~ 38.00°C (95.00°F ~ 100.40°F): ±0.05°C(0.09°F) 38.01°C ~ 45.00°C (100.42°F ~ 113.00°F): ±0.10°C(0.18°F)
Operating Condition	5 °C ~ 40 °C; Humidity 15% ~ 93%
Storage Condition	-25 °C ~ 70 °C; Humidity 10% ~ 95%
Resolution	0.01°C / 0.01°F
Measurement Unit	°C or °F Switchable
Power Source	CR2032
Battery Life	200 hours
Communication	Bluetooth 4.0
Bluetooth Distance	Maximum 10M
Dimension	40.5 (L) x 27.5 (W) x 11.7 (H) mm
Weight	18 g
Medical-grade Band	For upper arm circumferences: 12~20cm (4.7"~7.9")

SHIPPING INFO

Units per Carton	100
Gross Weight per Carton	9 KG
Carton size	40.3 (L) x 33 (W) x 32.1 (H) cm



MbH BM



TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288

www.taidoc.com