


GlucoX



Strip
automatic identification



GX1

Blood Glucose & β -Ketone MONITORING SYSTEM



NO CODING



MYLINK



LCD
BACKLIGHT



STRIP
EJECTION



EF TECHNOLOGY



BLUETOOTH

SPECIFICATIONS- Meter

Communication	Bluetooth or Strip port
Ketone Warning	Yes
Power Source	2 x AAA
Memory Capacity	1000 sets
Day Average	7, 14, 21, 28, 60, 90 days for blood glucose
Daily Alarm	4 daily alarms
Dimension	95 (L) x 64 (W) x 24 (H) mm
Weight	81g (without battery)
Operating Condition	10°C (50°F) ~ +40°C (104°F), 10% ~ 85% R.H.
Storage Condition	-20 °C ~ +60 °C (Meter); 2°C ~ 30°C (Strip)

SPECIFICATIONS- Strip

	Glucose	Ketone
Strip Type	Glucose	Ketone
Enzyme Type	GDH-FAD	HBD
Sample Size	0.5 μ L	1.0 μ L
Reaction Time	5 seconds	5 seconds
Measurement Range	20 ~ 650 mg/dL (1.1 ~ 36.1 mmol/L)	0.1 ~ 8.0 mmol/L
Hematocrit Range	15% ~ 65%	35% ~ 50%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at 100 mg/dL (5.55 mmol/L)	\leq 1 mmol/L, SD < 0.1 mM; > 1 mmol/L, CV < 7.5%
Accuracy	ISO 15197:2013 ; EN ISO 15197:2015	-
Package	Vial pack / Single foil pack	Vial pack / Single foil pack