

TD-4116


BLOOD GLUCOSE MONITORING SYSTEM

7 Sec, 0.7 μ L

Display Screen

FEATURES

- Strip indication light
- Large LCD with back light
- Before / After meal (AC/PC) presetable
- Ketone warning



SPECIFICATIONS

Enzyme Type	GOD
Sample Size	0.7 μ L
Reaction Time	7 seconds
Measurement Range	20 ~ 600mg/dL (1.1 ~ 33.3mmol/L)
Hematocrit Range	30% ~ 55%
Precision	SD < 5mg/dL (0.278mmol/L) at < 100mg/dL (5.55mmol/L); CV < 5% at 100mg/dL (5.55 mmol/L)
Accuracy	± 15 mg/dL if < 100mg/dL : $\pm 15\%$ if ≥ 100 mg/dL
Ketone Warning	Yes
Communication	Mini USB
Power Source	AAA x 1
Memory Capacity	450 sets
Day Average	7, 14, 21, 28, 60, 90 Days
Daily Alarm	4 daily alarms
Dimension	89.8 (L) x 52.4 (W) x 17.3 (H) mm
Weight	46.6 g (without battery)
Operating Condition	10°C (50°F) ~ 40°C (104°F)
Storage Condition	-20°C (-4°F) ~ +60°C (140°F) (Meter); 4°C (39.2°F) ~ 40°C (104°F) (Strip)
Units Per Carton	30
Carton Size	405 (L) x 390 (W) x 365 (H) mm
Gross Weight Per Carton	11.1 kg



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288