

Vet Ca²⁺

TD-5220

Plasma Calcium Monitoring System



FEATURES

- Rapid Serum Calcium Test
- Quick Plasma Separation

ACCESSORIES

Part no. 5220-CS001



Calcium test strip

SPECIFICATIONS

Sample Size	45 μ L \pm 5 μ L
Reaction Time	180 seconds
Measuring Range	0.5 to 3.0 mmol/L
Precision	\leq 1 mmol/L : SD < 0.12 mmol/L > 1 mmol/L : CV < 7.5%
Accuracy	\pm 25%
Communication	RS232
Power Source	AAA x 2
Memory Capacity	300 sets
Dimension	166.1 (L) x 34.9 (W) x 34.9 (H) mm
Weight	81 g (without battery)
Operating Condition	2°C ~ 35°C
Storage Condition	-20°C (-4°F) ~ +60°C (140°F) (Meter); 2°C (35.6°F) ~ 30°C (86°F) (Strip)

SHIPPING INFO

Units per Carton	50
Gross Weight per Carton	
Carton Size	38.8 (L) x 21.8 (W) x 37.5 (H) cm



TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.

www.taidoc.com