

CHOLESTEROL



GLUCOSE



KETONE



TRIPLE IN ONE

Use only with Procheck Witty Blood Glucose / β -Ketone / Total Cholesterol Test Strips

Procheck Smart Smartphone Dongle



Made for
iPhone | iPad | iPod



IOS Compatible Device (Requires operating system 14.0 or later)

1. iPad Air, iPad mini 3
 2. iPad (6th generation), iPod touch (6th generation)
 3. iPhone 5s / 6 / 6s Plus / 7 Plus / XR / XS / 11 / 12 Pro Max
- *Some ios phones might not be compatible



Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries and regions. App Store is a service mark of Apple Inc.

SPECIFICATIONS

Part No.	Procheck Smart
Communication	Lightning
Dimension	51 (L) X 24 (W) X 8.8 (H) mm
Weight	9 g
Operating Condition	+10°C (50°F) ~ +40°C (104°F)
Storage Condition	-20°C (-4°F) ~ +60°C (140°F)
Shipping INFO	
Units per Carton	44
Gross Weight per Carton	10.3 kg
Carton Size	467 (L) x 294 (W) x 415 (H) mm

SPECIFICATIONS / STRIP

GLUCOSE



Sample Size	0.8 μ L
Reaction Time	5 seconds
Measurement Range	10 ~ 700 mg/dL (0.56 ~ 38.89 mmol/L)
Hematocrit Range	0 ~ 70 %
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at \geq 100 mg/dL (5.55 mmol/L)
Package	Vial Pack (50 pcs x 1 vial)
Package Dimension	43 (L) x 40 (W) x 58 (H) mm
Shipping INFO	
Units per Carton	400
Gross Weight per Carton	15.9 kg
Carton Size	444 (L) x 393 (W) x 402 (H) mm

KETONE



Sample Size	1.0 μ L
Reaction Time	10 seconds
Measurement Range	0.1 ~ 8.0 mmol/L
Hematocrit Range	10 ~ 70 %
Precision	\leq 1 mmol/L, SD < 0.1 mM; > 1 mmol/L, CV < 7.5%
Package	Single Foil Pack (50 pcs)
Package Dimension	96 (L) x 10 (W) x 220 (H) mm
Shipping INFO	
Units per Carton	720
Gross Weight per Carton	NA
Carton Size	456 (L) x 378 (W) x 368 (H) mm

CHOLESTEROL



Sample Size	3.0 μ L
Reaction Time	\leq 60 seconds
Measurement Range	100 ~ 400 mg/dL (2.6 ~ 10.4 mmol/L)
Hematocrit Range	20 ~ 60%
Precision	CV < 7.5%
Package	Single Foil Pack (10 pcs)
Package Dimension	96 (L) x 5.5 (W) x 220 (H) mm
Shipping INFO	
Units per Carton	280
Gross Weight per Carton	NA
Carton Size	456 (L) x 378 (W) x 368 (H) mm