



CHOLESTEROL

GLUCOSE

KETONE

TRIPLE IN ONE

Use only with Procheck Witty Blood Glucose / β -Ketone / Total Cholesterol Test Strips

Procheck Witty

Smartphone Dongle



- USB-C Connector**
Connects directly to your Smartphone
- Pocket-Size**
No need for a bulky case
- ProCheck APP**
ProCheck app for personal health management
- Battery Free**
The power comes from the Smartphone

Android Compatible Device (Recommends operating system 12.0 or later)



1. Samsung Galaxy A70 / A8+ / S8 / S9 / S9+ / S10 / S10+ / Note 8 / Note 9 / Note 10
 2. Xiaomi Mi 9T / Mi 11 Lite / Redmi Note 7 / Redmi Note 8 Pro
 3. Google Pixel 3 / Pixel 4 / Pixel 5
 4. OPPO Reno Z 5. Realme XT 6. Nokia 8.1
- *Some android phones might not be compatible



Google Play and the Google Play logo are trademarks of Google LLC.

SPECIFICATIONS

| | |
|-------------------------|-----------------------------------|
| Part No. | Procheck Witty |
| Communication | USB-C |
| Dimension | 55 (L) X 40.29 (W) X 10.35 (H) mm |
| Weight | 11.8 g |
| Operating Condition | +10°C (50°F) ~ +40°C (104°F) |
| Storage Condition | -20°C (-4°F) ~ +60°C (140°F) |
| Shipping INFO | |
| Units per Carton | 44 |
| Gross Weight per Carton | 10.3 kg |
| Carton Size | 467 (L) x 294 (W) x 415 (H) mm |

SPECIFICATIONS / STRIP

GLUCOSE



| | |
|-------------------------|--|
| Sample Size | 0.8 μ L |
| Reaction Time | 5 seconds |
| Measurement Range | 10 ~ 700 mg/dL (0.56 ~ 38.89 mmol/L) |
| Hematocrit Range | 0 ~ 70 % |
| Precision | SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at \geq 100 mg/dL (5.55 mmol/L) |
| Package | Vial Pack (50 pcs x 1 vial) |
| Package Dimension | 43 (L) x 40 (W) x 58 (H) mm |
| Shipping INFO | |
| Units per Carton | 400 |
| Gross Weight per Carton | 15.9 kg |
| Carton Size | 444 (L) x 393 (W) x 402 (H) mm |

KETONE



| | |
|-------------------------|---|
| Sample Size | 1.0 μ L |
| Reaction Time | 10 seconds |
| Measurement Range | 0.1 ~ 8.0 mmol/L |
| Hematocrit Range | 10 ~ 70 % |
| Precision | \leq 1 mmol/L, SD < 0.1 mM; > 1 mmol/L, CV < 7.5% |
| Package | Single Foil Pack (50 pcs) |
| Package Dimension | 96 (L) x 10 (W) x 220 (H) mm |
| Shipping INFO | |
| Units per Carton | 720 |
| Gross Weight per Carton | NA |
| Carton Size | 456 (L) x 378 (W) x 368 (H) mm |

CHOLESTEROL



| | |
|-------------------------|-------------------------------------|
| Sample Size | 3.0 μ L |
| Reaction Time | \leq 60 seconds |
| Measurement Range | 100 ~ 400 mg/dL (2.6 ~ 10.4 mmol/L) |
| Hematocrit Range | 20 ~ 60% |
| Precision | CV < 7.5% |
| Package | Single Foil Pack (10 pcs) |
| Package Dimension | 96 (L) x 5.5 (W) x 220 (H) mm |
| Shipping INFO | |
| Units per Carton | 280 |
| Gross Weight per Carton | NA |
| Carton Size | 456 (L) x 378 (W) x 368 (H) mm |