

TD-4277

BLOOD GLUCOSE MONITORING SYSTEM



NEW

3 Button Design


Easy to review the test results

Indicator light for target range

Individual target range setting



Apple and the Apple logo are trademarks of Apple Inc.



Google Play and the Google Play logo are trademarks of Google LLC.



SPECIFICATIONS

Part No.	TD-4277C
Enzyme Type	GDH-FAD
Sample Size	0.5 μ L
Reaction Time	5 seconds
Measurement Range	20 ~ 600 mg/dL (1.1 ~ 33.3 mmol/L)
Hematocrit Range	20% - 65%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at 100 mg/dL (5.55 mmol/L)
Accuracy	\pm 15 mg/dL if < 100 mg/dL; \pm 15% if \geq 100 mg/dL
Ketone Warning	Yes
Communication	Bluetooth or USB Type C
Power Source	AAA x 2
Memory Capacity	1000 sets
Day Average	7, 14, 21, 28, 60, 90 days
Daily Alarm	4 daily alarms
Dimension	96 (L) x 61 (W) x 26 (H)mm
Weight	67.2 g (without battery)
Operating Condition	10°C (50°F) ~ 40°C (104°F), below 85% R.H.
Storage Condition	-20°C (-4°F) ~ +60°C (140°F) (Meter) 2°C (35.6°F) ~ 32°C (89.6°F) (Strip)

SHIPPING INFO

Units per Carton	40
Gross Weight per Carton	12.8 kg
Carton Size	49.8 (L) x 35.8 (W) x 36 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.