

TD-4125

BLOOD GLUCOSE MONITORING

7 Sec, 0.7 μ L

FEATURES

- Adjustable Hi/Lo indicator
- LCD Backlight
- Ketone warning



 **Bluetooth**[®]



Download on the
App Store

Apple and the Apple logo
are trademarks of Apple Inc.




GET IT ON
Google Play

Google Play and the Google
Play logo are trademarks of
Google LLC.



SPECIFICATIONS

Enzyme Type	GOD
Sample Size	0.7 μ L
Reaction Time	7 seconds
Measurement Range	20 ~ 600 mg/dL (1.1 ~ 33.3 mmol/L)
Hematocrit Range	30% - 55%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at 100 mg/dL (5.55 mmol/L)
Accuracy	\pm 15 mg/dL if < 100 mg/dL; \pm 15% if \geq 100 mg/dL
Ketone Warning	Yes
Communication	RS232-4 Poles or Bluetooth
Power Source	AAA x 2
Memory Capacity	450 sets
Day Average	7, 14, 21, 28, 60, 90 days
Daily Alarm	4 daily alarms
Dimension	97 (L) x 52 (W) x 16 (H) mm
Weight	47.8 g (without battery)
Operating Condition	10°C (50°F) ~ 40°C (104°F), below 85% R.H.
Storage Condition	-20°C (-4°F) ~ 60°C (140°F) (Meter); 4°C (39.2°F) ~ 40°C (104°F) (Strip)
Units per Carton	48
Carton Size	45.6 (L) x 37.8 (W) x 36.8 (H) cm
Gross Weight per Carton	12.8 kg



www.taidoc.com

TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.