

TD-5212

Digital Single Use Pregnancy Test



FEATURES

- Easy to use
- Result in 5 minutes
- Clinical Accuracy > 99%



SPECIFICATIONS

Model Number	TD-5212A
hCG Detection Level	25 mIU
Device Dimension	141 x 26.5 x 16.2 mm
Weight	24 g
Operational Condition	In Room Temperature
Storage Condition	2 °C ~ 30 °C with RH lower than 85%
Units per carton:	nn Units
Net weight per carton:	nn kg
Gross weight per carton:	nn kg



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel : +886-2-6625-8188 Fax : +886-2-6625-0288