

# TD-5302

## TRADITIONAL PREGNANCY TEST



### FEATURES

- 99% accurate at detecting typical hormone levels
- Reaction time less than 5 minutes
- One step easy operation



### SPECIFICATIONS

hCG Detection Level	25 mIU
Device Dimension	131 x 23.2 x 8.3 mm
Weight	11.5 g
Operational Condition	In Room Temperature
Storage Condition	2°C (35.6°F) ~ 30°C (86°F), with RH lower than 85%
Unit per Carton	45
Gross Weight per Carton	3.25 kg
Carton Size	L31.8 x W29.7 x H21.5 cm



TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,  
New Taipei City 24888, Taiwan  
Tel: +886-2-6625-8188 Fax: +886-2-6625-0288