



TD-2555

WEIGHT SCALE



Body mass index function



Age, height, and gender for more precise algorithm



Large LCD with backlight



Automatically power-off function



iOS



Android



Healthy Check



iOS



Android



Healthy Check Pro

SPECIFICATIONS

BMI (Body Mass Index)	Yes
Max. capacity	250 kg
Accuracy	4.0 kg ~ 60.0 kg $\pm$ 0.3 kg, 60.1 kg ~ 250.0 kg $\pm$ 0.5%
Graduation	0.1 kg
Age Range	10 to 80
Height Range	100 to 220 cm
User	5
Communication	Bluetooth smart
Memory	135 sets
Backlight	LCD backlight
Power Source	4 x AAA
Dimension	340(L) x 400 (W) x 26 (H) mm

SHIPPING INFO

Units per Carton	5
Gross Weight per Carton	18 kg
Carton Size	49.3 (L) x 29.3 (W) x 45.5 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.