

# URiGHT

## EZTemp Infrared Thermometer



Infrared Technology




Magnetic Probe Cap  
Ear/Forehead Switchable



Fever Warning



1 sec Scan



Multi Measurement  
Modes

### SPECIFICATIONS



Accuracy	$\pm 0.2^{\circ}\text{C}$ ( $\pm 0.4^{\circ}\text{F}$ ) for the range of 35 - 42 ( $95^{\circ}\text{F}$ - $107.6^{\circ}\text{F}$ ) $\pm 0.3^{\circ}\text{C}$ ( $\pm 0.5^{\circ}\text{F}$ ) for the range of $< 35$ ( $95^{\circ}\text{F}$ ) or $> 42$ ( $107.6^{\circ}\text{F}$ )
Measurement range	$32.0^{\circ}\text{C}$ - $42.9^{\circ}\text{C}$ ( $89.6^{\circ}\text{F}$ - $109.2^{\circ}\text{F}$ )
Operation conditions	$10^{\circ}\text{C}$ - $40^{\circ}\text{C}$ ( $50^{\circ}\text{F}$ - $104^{\circ}\text{F}$ ) 15 - 95% non-condensing
Resolution	$0.1^{\circ}\text{C}/0.1^{\circ}\text{F}$
Display	LCD with 4-color backlit, White/Green/Orange/Red
Backlit interpretation	Power on: White $32^{\circ}\text{C} \leq T \leq 37.5^{\circ}\text{C}$ $37.6^{\circ}\text{C} \leq T \leq 38^{\circ}\text{C}$ $38.1^{\circ}\text{C} \leq T \leq 42.9^{\circ}\text{C}$
Communication	Bluetooth 4.0
Measurement mode	Forehead, Ear and Object
Measurement unit	$^{\circ}\text{C}$ or $^{\circ}\text{F}$ switchable
Magnetic Probe Cap	Yes
Power Source	AAA x 2
Battery life	At least 3,000 measurements
Dimension	140mm x 35mm x 40mm
Weight	60g (w/o batteries)
Packing	100 PCS/CTN/11KGs, 50x32x27cm