

L1 Lactate Monitoring System

FEATURES

- XPER strip technology
- Expiry date reminder



Apple and the Apple logo are trademarks of Apple Inc.



Google Play and the Google Play logo are trademarks of Google LLC.



SPECIFICATIONS

Part No.	TD-4261C
Reaction time	5 seconds
Measurement Range	0.3 - 22 mmol/L
Hematocrit Range	10% - 65%
Sample Size	0.8 μ L
Precision	≤ 3 mmol/L, $SD < 0.3$ mM : > 3 mmol/L, $CV < 7.5\%$
Package	Single foil pack
Communication	Bluetooth or USB Type-C
Power Source	AAA x 2
Memory Capacity	1000 sets
Dimension	93(L) x 57(W) x 18(H) mm
Weight	50.2g (without battery)
Operating Condition	10°C (50°F) ~ 40°C (104°F), 10~85% RH
Storage Condition	-20°C to 60°C (-4°F to 140°F) (Meter); 2°C to 30°C (35.6°F to 86°F) (Strip)

Shipping Info

Units Per Carton	42
Gross Weight Per Carton	13 kg
Carton Size	45.6 (L) x 37.8 (W) x 36.8 (H) cm

Accessories

ACC-006	Type-A to Type-C USB cable
---------	----------------------------