

TD-1090

Temperature Chamber

FEATURES

- Accuracy verification for ear and forehead thermometers
- Adjustable target temperature
- Data export capability



Thermometer Fixter

Communication Cable

Part no.	ACC-1090-001	ACC-1090-002	ACC-1090-003	ACC-1090-004	ACC-1090-005
Picture	 TD-1107	 TD-1261	 EZTemp	 TD-1241	 TD-1242

SPECIFICATIONS

Display Resolution	0.1°C/0.1°F
Power & Fuse(s)	100 ~ 240VAC (50/60Hz) / 3A
Memory	1000
Operating Condition	16°C to 32°C (60.8°F to 89.6°F), 85% RH or less
Temperature Range (0.98 emissivity)	37.0°C ~ 41.0°C (98.6°F ~ 105.8°F)
Acceptable Recalibration Range	±0.5°C (±0.9°F)
Accuracy	±0.1°C
Stability	±0.03°C at 37°C
Measurement Units	°C (default) or °F
Warm-up Time at 25°C Ambient	30 minutes
Nominal Emissivity	0.98 emissivity
Computer Interface	USB HID
Dimension	300mm x 350mm x 160mm

SHIPPING INFO

Units per Carton	1
Gross Weight per Carton	6.93 kg
Carton size	41.3 (L) x 22.3 (W) x 37.6 (H) cm



MT-1099

CLEVER Reference Meter



Display resolution	0.01°C
Measurement range	32°C-43°C
Accuracy	±0.05°C
Operating condition	16°C-32°C, 30%-85%RH
Measurement unit	°C
Traceability	Trace to National Measurement Laboratory R.O.C
Calibration interval	1 year
Battery life	>1,000 measurements

Fixter

Part no.: ACC-1090-002

