

URiGHT

ADVANCE WEIGHT SCALE

For Baby, Pet and Adult



BABY GROWTH TRACKING

PET HEALTH MANAGEMENT

ADULT DAILY MONITORING

- Automatic power on/off
- Low battery / Overload indication
- Four High precision strain gauge sensors
- Easy to assemble



SPECIFICATIONS

Part No.	TD-2591
Max. Capacity	180 kg
Capacity	50g~50kg (d=10g), 50~150kg (d=50g), 150kg~180kg (d=100g)
Unit	kg, lb, st:lb
Backlight	LCD backlight
Power Source	3 x AAA
Dimensions	300(L) x 300(W) x 25(H) mm
Dimensions with tray	635(L) x 307.5(W) x 55 (H) mm

SHIPPING INFO

Units per Carton	5
Gross Weight per Carton	13.5 kg
Carton Size	50.5 (L) x 38.5 (W) x 36 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.