


TD-2555




4G Cellular Weight Scale

4G connectivity via SIM card, with eSIM available as an option



Telemedicine without time difference!



4G MODULE

Part. No.	Description	Specification
EXS62-W	LTE 4G module for American Requirement	Frequency bands: 700 MHz (Bd12), 850 MHz (Bd5), AWS (Bd4), 1900 MHz (Bd2)



TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.



Body mass index function



Age, height, and gender for more precise algorithm



Large LCD with backlight



Automatically power-off function

SPECIFICATIONS

BMI (Body Mass Index)	Yes
Max. capacity	250 kg
Accuracy	4.0 kg ~ 60.0 kg \pm 0.3 kg, 60.1 kg ~ 250.0 kg \pm 0.5%
Graduation	0.1 kg
Age Range	10 to 80
Height Range	100 to 220 cm
User	5
Communication	4G
Memory	135 sets
Backlight	LCD backlight
Power Source	4 x AAA
Dimension	340(L) x 400 (W) x 26 (H) mm

SHIPPING INFO

Units per Carton	5
Gross Weight per Carton	18 kg
Carton Size	49.3 (L) x 29.3 (W) x 45.5 (H) cm



TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation.
All rights reserved.