

# URiGHT

## EZTemp Infrared Thermometer



Infrared Technology




Magnetic Probe Cap  
Ear/Forehead Switchable



Fever Warning



1 sec Scan



Multi Measurement  
Modes

### SPECIFICATIONS



Part No.	AOJ-20A
Accuracy	±0.2°C (±0.4°F) for the range of 35 - 42 (95°F - 107.6°F) ±0.3°C (±0.5°F) for the range of < 35 (95°F) or > 42 (107.6°F)
Measurement range	32.0°C - 42.9°C (89.6°F - 109.2°F)
Operation conditions	10°C - 40°C (50°F - 104°F) 15 - 95% non-condensing
Resolution	0.1°C/0.1°F
Display	LCD with 4-color backlit, White/Green/Orange/Red
Backlit interpretation	Power on: White <div style="display: flex; justify-content: space-around; align-items: center;"> <span style="background-color: #90EE90; padding: 2px;">32.0°C ≤ T ≤ 37.5°C</span> <span style="background-color: #FFA500; padding: 2px;">37.6°C ≤ T ≤ 38°C</span> <span style="background-color: #FF0000; padding: 2px;">38.1°C ≤ T ≤ 42.9°C</span> </div>
Communication	Bluetooth 4.0
Measurement mode	Forehead, Ear and Surface
Measurement unit	°C or °F switchable
Magnetic Probe Cap	Yes
Power Source	AAA x 2
Battery life	At least 3,000 measurements
Dimension	140mm x 35mm x 40mm
Weight	60g (w/o batteries)
Packing	100 PCS/CTN/11KGs, 50x32x27cm