


# URIGHT

## ADVANCE Arm Blood Pressure Monitor



Bluetooth



Irregular Heartbeat (IHB) Detection



Cuff Wearing Detection



Keep Still indicator



iOS

Android



Healthy Check

### SPECIFICATIONS

Part No.	AOJ-30B
Accuracy of Pressure	±3 mmHg
Accuracy of Heart Rate	±5 bpm
Measurement Range	Systolic : 57 ~ 255 mmHg / Diastolic : 25 ~ 195 mmHg Pulse rate : 40~199 beat per minute
Memory Capacity	198 memory with date and time
Communication	Bluetooth 4.0
Power Source	4 x AAA DC6V or DC5V type C cable(Optional)
Dimension	127 (L) x 93 (W) x 73 (H) mm
Weight	262g (without battery)
Operating Condition	5°C~ 40°C ; Humidity 15% ~ 95%
Battery Life	At least 200 measurements
Power Saving	60 seconds (with any operation)
Units per Carton	24
Gross Weight per Carton	
Carton Size	448 (L) x 370 (W) x 348 (H) mm

### ACCESSORY

Cuff Size	Arm Circumference
Regular	22-42 cm



www.taidoc.com

TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.