

TD-3128

BLOOD PRESSURE MONITORING SYSTEM



Data Transmission

4 Users Selection

FEATURES



• Average 3 measurements



• Large screen 121 x 63mm



• IHB Detection



• Low-noise



iOS

Android



Healthy Check

SPECIFICATIONS

Accuracy of Pressure	± 3 mmHg or $\pm 2\%$ of reading
Accuracy of Heart Rate	$\pm 4\%$ of reading
Measurement Range	Systolic : 60 ~ 255 mmHg / Diastolic : 30 ~ 195 mmHg Pulse rate : 40~199 / minute
Memory Capacity	400 sets
Data Output	Bluetooth 4.0 (TD-3128B) / RS232(TD-3128C)
Multi-user	4
Power Source	4 x AA or 6V DC Adaptor (Optional)
Dimension	141 (L) x 121 (W) x 72 (H) mm
Weight	350 g (without battery)
Operating Condition	5 °C ~ 40 °C ; Humidity 15% ~ 93%
Battery Life	At least 200 measurements
Power Saving	180 seconds

ACCESSORY

Part No.	Cuff Size	Arm Circumference
RV01	Small	15~24 cm
RV02	Regular	24~35 cm
RV03	Large	32~44 cm
RV04	Wide Range	24~43 cm



SHIPPING INFO

Units per Carton	12
Gross Weight per Carton	12 kg
Carton Size	43.3 (L) x 37.8 (W) x 36.5 (H) cm



www.taidoc.com

TaiDoc Technology Corp.

Copyright TaiDoc Technology Corporation. All rights reserved.